

Bronocide™ STERILISING SOLUTION

LASTING
BIOCIDAL
ACTION

Effective sterilisation of all instruments & surfaces.

READY FOR USE - DO NOT DILUTE

BACTERIACIDAL • FUNGICIDAL • SPORICIDAL • VIRUCIDAL



Available in: 1 Litre, 5 Litre and 25 Litre

USE:

Effective STERILISATION of all surgical stainless steel instruments where autoclaving is not possible.
SAFE STERILISATION of Endoscopes and related equipment. STERILISATION of all medical plastic devices and tubes.
Sterile facilities and clean rooms in the healthcare and pharmaceutical industries. Food and dairy processing plants.
In-vitro fertilisation facilities. Laboratories in general. Warehouse and storage facilities for perishable goods.
Plant and tissue culture facilities. Indoor and outdoor commercial horticultural / agricultural enterprises.

NOTE:

BRONOCIDE Sterilising solution is a highly effective general sterilising agent with the same efficacy in airspace fogging applications in all sensitive areas where micro-biological contamination poses a problem.
Non-Lacramatic (No burning of eyes)

INSTRUCTIONS FOR USE:

1. Pour solution into suitable container with lid.
2. Immerse instruments/items to be sterilised in the solution.
3. Completely immerse for required time (5-15min).
4. DO NOT RINSE - when removing items for use. Shake off excess solution and use directly.
5. Stainless steel (316) items can be left fully immersed for extended periods.
6. Non-stainless steel items can be immersed for up to 4 hours without experiencing corrosion issues.

Active ingredient 2%

Efficacy of the disinfectant could be compromised if surfaces are soiled.

CONTAINS

Dichlorobenzyl alcohol, Bromonitropropane diol, Isopropyl alcohol, Benzalkonium chloride.

CAUTION

Keep out of reach of children. In case of accidental eye contact, rinse with plenty of water.

STORE BELOW 25°C

NRCS Reg. No. Act5GNR529/290979/090/0927

FIRST AID MEASURES

Eyes: Wash immediately with copious amounts of water and obtain medical attention.

Skin: Prolonged contact may cause irritation with sensitive skin types. Avoid prolonged skin contact and wash off with water if irritation occurs.

Ingestion: Rinse mouth and drink 2 glasses of water or milk. If a large amount is ingested and symptoms develop, seek medical attention. Do not induce vomiting.

PERSONAL PROTECTIVE EQUIPMENT:

Eyes: When the product is used as directed, personal protective equipment is not required. However, the product is an eye irritant, and as such it should be handled in such a way as to avoid eye contact.

Hand/Skin: Prolonged skin contact should be avoided, and if an application dictates otherwise, the use of protective gloves or clothing is recommended.

Respiration: No special protection is required.